Exclusive Distribution by Plasma Trading Joint-Stock Company (PLD)

HEAD OFFICE

Lot 06, B row, TT19, Van Quan Urban Area, Phuc La Ward, Ha Dong Dist., Ha Noi

Mobile: 024.6294.7733

Hotline: 0901702626

Email: info@plasmamed.com.vn

www.plasmamed.com.vn

HCM REPRESENTATIVE OFFICE

No.6, 5A Cu Xa Binh Thoi St, Ward

Mobile: 028.3962.2379

Hotline: 08.8899.2222

Email: info.hcm@plasmamed.com.vn

Website: www.plasmamed.com.vn

DA NANG REPRESENTATIVE OFFICE

No.68 Duong Ba Trac, Hoa Cuong Nam Ward, Hai Chau Dist, Da Nang

Mobile: 098.255.8114

Email: info.dn@plasmamed.com.vn

Website: www.plasmamed.com.vn

VIETNAM PLASMA TECHNOLOGY JOINT STOCK COMPANY

CATALOGUE

THE COLD PLASMA

TECHNOLOGICAL PIONEERS IN THE WORLD

Company Information

MANUFACTURER: Vietnam Plasma Technology Joint-Stock Company (PLT)

Office Address: A4 - TT6, **Van Quan Urban Area, Phuc La Ward, Ha Dong Dist., Hanoi**

Telephone: 024 6294 7733

Email: info@plasma.vn

Website: www.plasma.vn

EXCLUSIVE DISTRIBUTOR: Plasma Trading Joint-Stock Company (PLD)

Website: www.plasmamed.com.vn

www.plasma.vn

Telephone: 024.6294.7733

Hotline: 090.170.2626

Email: info@plasmamed.com.vn;

info@plasma.vn

Address:

- **Head Office:** Lot 06, B row, TT19, Van Quan Urban Area, Phuc La Ward, Ha Dong Dist., Hanoi
- **HCM Representative:** No.6, 5A Cu Xa Binh Thoi St., Ward 08, Dist 11., Ho Chi Minh
- **Da Nang Representative:** No.68 Duong Ba Trac, Hoa Cuong Nam Ward, Hai Chau Dist, Da Nang

| Vietnam Plasma Technology Joint-Stock Company

INTRODUCTION

PlasmaMED® is a cold plasma beam generator successfully invented and manufactured by Plasma Vietnam Technology Corporation (PLT). It effectively kills bacteria, including antibiotics resistant bacteria, and improves wound healing by stimulating tissue regeneration, skin cell proliferation and neovascular proliferation without causing any side effects. PlasmaMED® not only reduces time and cost for treating slow-healing wounds but also is an effective tool for treating open wounds in general and a number of skin diseases. For post-operative wounds, treatment with


Antibiotics resistant bacteria makes wound treatment more costly and prolonged. According to the World Health Organization, antibiotic resistance has become an alarming threat, especially in developing countries.

Slow healing of wounds may be caused by many other factors, including but not limited to aging, chronic conditions, microvascular injury and immunodeficiency.

Slow healing of wounds may result in complications and thus higher risks of mortality, imposing treatment burdens on patients, their families and hospitals.

Current versions of PlasmaMED® machine

Model: PlasmaMED-GAP


■ Shell size:

Height	Width	Depth	Material
960 mm	405 mm	530 mm	ABS Plastic

■ The size of Machine - Base

Height	Width	Depth	Material
50 mm	490 mm	490 mm	ABS Plastic

Model: PlasmaMED-TM

■ Shell size:

Height	Width	Depth	Material
910 mm	255 mm	480 mm	ABS Plastic

■ The size of Machine - Base

Height	Width	Depth	Material
50 mm	355 mm	505 mm	ABS Plastic

Some Clinical Results

1. Treatment of chronic heel ulcer due to osteitis, antibiotic resistance:

Case description:

Female, aged 42, chronic heel ulcer due to osteitis, infected with bacteria resistant to all antibiotics, wound dehiscence, remaining unhealed after one month of normal wound dressing changes. The results for clinical examination by Thu Duc District Hospital

Course of treatment:


Patient was treated with radiation of cold plasma beam from PlasmaMED®. Radiated the whole wound and its adjacent areas. Started from the deepest spot, then gradually moved to the outer parts at 5 mm/s speed. Normal wound dressing change on daily basis, no additional medication.

Treatment outcome:

After 6 sessions of radiation, wound healed, no longer any sign of infection, patient checked out hospital for family-based aftercare.

Doctor's comment:

Very good outcome. The wound would not have healed if not treated with cold plasma.

chronic heel ulcer due to osteitis

Cutting operations

1st treatment


2nd treatment

3rd treatment

4th treatment

5th treatment

6th treatment

After 8 day treatment

Follow up 6 months

2. Treatment of eczema with PlasmaMED®:

Case description:

Female, aged 33, lesion caused itching and papules on the front of inner side of right lower leg. First lesion appeared 7 years ago, patient experienced multiple periods of itching and papules ever since. Diagnosed with eczema and treated at the National Hospital of Dermatology and Venereology. Patient was treated with corticoid cream, anti-histamine, and anti-biotics during period of super-infection. Lesions went away, but came back after some time with increased frequency and in larger area.

Course of treatment:

Patient was treated with cold plasma from PlasmaMED® with average radiation duration of 20s/cm², twice a day in the first week and then once a day in the second week. Treatment finished after 2 weeks.

Treatment outcome:

No more itching in lesions around the primary lesion and skin surface returned to normal color. Surface of the primary lesion returned to flat state, no more itching. Follow-up examination after 8 months showed no re-appearance of lesion.

Doctor's comment:

Very good results, the wound would not have healed if not treated with PlasmaMED®.

3. Treatment of pressure ulcer:

Case description:

Female, aged 87, had stroke and then had to stay in bed as a consequence. Trochanter major bed sore acquired after one month in bed. Trochanter minor bed sore acquired early the third month in bed. Excess discharge of fluid, blood and greenish pus and no sign of healing. Bed sores got larger and larger, causing pain, poor eating and sleeping, and significant weight loss to patient. The results for clinical examination by **National Hospital of Endocrinology**.

Course of treatment:

Patient treated with cold plasma from PlasmaMED® with average radiation duration of 20s/cm², frequency once a day. Treatment completed after 30 days.

Treatment outcome:

The ulcers improved significantly day by day. After 5 days of radiation, ulcer surface got dry, discharge stopped, ulcers began to contract, epithelialization started. After 30 days of radiation therapy, ulcers closed completely and epithelialization completed.

Doctor's comment:

Have very good results, **the wound was healed positively**

(01/07/17 - Mấu trái) - The left

(01/07/17 - Mấu phải) - The right

(31/07/17 - Mấu trái) - The left

(31/07/17 - Mấu phải) - The right

(03/07/16)

(10/07/16)

4. Treatment of ulcer caused by wound infection:

Case description:

Male, aged 36, suffered a wood saw injury in his thumb and forefinger of left hand. Wound got infected and ulcer got increasingly larger. Already treated with antibiotics for one month but excess exudate continued, ulcer area continued expanding. The results for clinical examination by Hue Central Hospital

Course of treatment:

Patient was treated with cold plasma from PlasmaMED® with average radiation duration of 20s/cm², frequency once a day. Treatment finished after 7 days.

Treatment outcome:

Wound improved significantly day by day. After 3 days of radiation, ulcer surface got dry, wound exudate stopped, ulcer surface showed light pink color. Epitheliazation completed after 7 days of radiation.

Doctor's comment:

Have very good results by plasma, the wound was healed positively, and appeared granulation.

5. Toe necrosis is resulted from complications of diabetes type 2:

Case description:

Male patient, 72 years old. In April, 2017, his I-II toes (P) are ulcerated, and have no sign of healing, leading to necrosis.

Patient has just discharged from Nguyen Tri Phuong Hospital – Endocrinology Department on July 3rd, 2017 with the state: His I-II toes are in dry necrosis, and he is living with diabetes type 2.

Then, patient gets the more and more painful at the I-II toes. His II toe releases pus which is full of stench. He forces to be hospitalized again on July 17th, 2017.

Course of treatment:

Patient is radiated cold plasma by PlasmaMED®. The whole of wound and its surrounding areas are rayed, from the deepest point then to the surrounding ones with the speed 5 mm/s. The gauze pads at the wound should be changed daily.

Use prescription: Natri clorid 0.9% 500ml, Paracetamol 37.5mg, Tramadol 325mg, Clopidogrel 75mg, Rosuvastatin 10mg, Telmisartan 40mg, Ceftriaxon 1g, Natri clorid

Treatment outcome:

After 5 times of treatment, the wound is healed. The patient is discharged from hospital and looked after at his home.

Doctor's comment:

Have very good results, because the wound had no sign of healing before.

Photos of treatment:

On July 19th, 2017: Getting the status of disease

Wash the wound and cut-off, filter

Radiate by the first time (2 minutes)

Radiate by the second time (2 minutes)

Radiate by the third time (2 minutes)

Radiate the fifth time (1 minutes)

– To be discharged from hospital

Photos of treatment:

Radiate by the first time

Radiate by the 5th time

Radiate by the 7th time

Radiate by the 9th time

Radiate by the 11th time

6. The wound is caused by bash:

Case description:

Male patient, 53 years old. In April , 2017, his wrist was injured (T), infected , appeared pus, swelled and painful. After the wound had no sign of healing for a month, he forced to be hospitalized .

He was treated at Cho Ray Hospital, then moved to Ho Chi Minh Orthopedics and Rehabilitation Hospital to receive the treatment .

Course of treatment:

Patient was radiated cold plasma by PlasmaMED. The whole of wound and its surrounding areas were rayed, from the deepest points then to the surrounding ones with the speed 5mm/s. The gauze pads at the wound should be changed daily.

Treatment outcome:

After 11 times of treatment, the wound is healed. The patient was discharged from hospital and looked after at his home.

Doctor's comment:

After being treated by plasma, the wound was healed positively, and appeared granulation. Following up examination after being discharged from hospital for a week, the wound was healed fully, and his wrist could move well.

Use of PlasmaMED[®] in Hospitals in Vietnam

Northern Vietnam:

1. National Hospital of Endocrinology
2. National Children's Hospital
3. National Hospital of Obstetrics and Gynecology
4. Military Hospital 103
5. National Hospital of Dermatology and Venereology
6. National Geriatric Hospital
7. Viet Duc University Hospital

Central Vietnam:

1. Hue Central Hospital
2. Thanh Hoa General Hospital
3. Thanh Hoa Children's Hospital
4. Hospital of Ngoc Lac (Thanh Hoa Province)

Southern Vietnam:

1. University Medical Center (Ho Chi Minh City HCMC)
2. Cho Ray Hospital (HCMC)
3. Gia Dinh Hospital (HCMC)
4. Trung Vuong Hospital (HCMC)

6. HCMC Hospital for Traumatology and Orthopaedics
7. Orthopaedics and Rehabilitation Hospital (HCMC)
8. Cu Chi General Hospital (HCMC)
9. Cu Chi District Hospital (HCMC)
10. Hospital for Nursing, Rehabilitation and Treatment of Occupational Diseases (HCMC)
11. People's Hospital 115 (HCMC)
12. General Thong Nhat Dong Nai Hospital
13. Dong Nai General Hospital (Dong Nai Province)
14. Ba Ria Hospital (Ba Ria-Vung Tau Province)
15. Nguyen Dinh Chieu General Hospital (Ben Tre)

CERTIFICATE OF LEGALITY

Certificate Registration For Circulation Of Medical Devices Manufacturing In Viet Nam

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CERTIFICATE

REGISTRATION FOR CIRCULATION OF
MEDICAL DEVICES MANUFACTURING IN VIETNAM

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CERTIFICATE
REGISTRATION FOR CIRCULATION OF
MEDICAL DEVICES MANUFACTURING IN VIETNAM

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- Pursuant to Decree No. 63/2012/ND-CP dated August 31st, 2012 issued by Government stipulating the functions, tasks, authority and organizational structure of the Ministry of Health;

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- Pursuant to Law on Quality of products and goods dated November 21st, 2007;

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- Pursuant to Circular No. 07/2002/TT-BYT dated May 30, 2002 of the Ministry of Health on guiding for circulation registration of medical devices.

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- Having examination of documentation and application letter for circulation of medical device submitted by the applicant.

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- Certificate of Science and Technology

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VN 14627

Moving Forward

Our research and innovation efforts have received support from Vietnam's Ministry of Health, Ministry of Science and Technology, and Hanoi Department of Science and Technology. Our aspiration is to become a regional leader in plasma technology with plasma application in various fields.

Building on the success of PlasmaMED®, we are developing cold plasma beam generators customized for treatment in different areas of specialization, such as pediatrics, dentistry, and gynecology.

In addition, we are undertaking research in application of cold plasma in treating contaminated food or chemical residues in food – a serious problem to the society these days. For a longer time horizon, we plan to explore application of plasma technology in agriculture, surface treatment industry, environmental technology and water treatment.

Support of and feedback from partners, scientists and regulators are an important source of encouragement to us. Our sincere thanks go to all of them!