



TANAPHAR®

SUPPLIER PORTFOLIO

TANAPHAR JOINT STOCK COMPANY

*Leading Medical Device & Healthcare
Manufacturer in Vietnam*

CONTACT INFORMATION

- **Address:** Thanh Oai Industrial Cluster, Binh Minh, Hanoi, Vietnam.
- **Primary Contact:** Alex Dao (Ms.)
- **Phone/WhatsApp/Viber/Telegram:** (+84) 968 174 279
- **Email:** Anhdao@tanaphar.com
- **Website:** Tanaphar.com

COVER LETTER

Dear Sir/Madam,

On behalf of **Tanaphar Joint Stock Company**, I would like to extend our warmest greetings and best wishes for your continued success.

With 20 years of dedication to the healthcare industry, Tanaphar stands as a leading manufacturer in Vietnam, specializing in high-quality medical consumables. Our state-of-the-art 15,000 m² facility is equipped with advanced automated technology from Japan and Germany, operating under strict ISO 13485:2016 and GMP standards.

We are proud to have successfully exported our products to many global markets. Our comprehensive product portfolio includes:

- Respiratory & Protection: Medical masks (3-layer, 4-layer, N95) and Herbal Cough Lozenges (effective for throat relief).
- Injection & Surgical: CE-certified syringes, needles, and sterile gauze.
- Wound Care: Medical tapes and UGOTANA adhesive bandages.
- In Vitro Diagnostics (IVD): Rapid Pregnancy Test strips and Ovulation Test kits.
- Personal & Pediatric Care: Fever Cooling Patches, Dental Floss Picks, safety nasal aspirators, and sterile Tongue Depressors/Cleaners.

We are looking for reputable partners to distribute our high-quality, compliant healthcare solutions. We also offer flexible OEM/ODM services tailored to your brand's needs.

Please find our Company Portfolio attached. I look forward to discussing how we can add value to your supply chain.

Sincerely,


Đào Thị Thùy Anh

Supply Chain & Import-Export Director




COMPANY OVERVIEW





- Name: Tanaphar Joint Stock Company
- Established: 2006
- Industry: Manufacturing and Distribution of Medical Disposables, Healthcare, and Personal Care Products.



Registration :
0101988259

**COMPANY
PROFILE**

COMPANY OVERVIEW

I. Production Capacity

Founded in 2007, Tanaphar is a leading manufacturer of medical consumables and healthcare products in Vietnam. Operating a 15,000 m² state-of-the-art facility with three independent specialized workshops, we are equipped with advanced automated production lines imported from Germany and Japan. Our entire manufacturing process strictly adheres to ISO 13485:2016 and GMP standards, ensuring that every product meets the highest international quality and safety requirements. With nearly two decades of expertise and a commitment to excellence, we currently



- **Workforce:** Over 200 skilled workers.
- **High Production Capacity:** 25,000,000 units per month, with the flexibility to scale up production immediately to meet market surges.
- **Automated Packaging Lines:** Utilizing high-speed Japanese packaging technology integrated with an end-to-end quality control system. This ensures superior hygiene, high productivity, and minimized manual

II. Financial Capacity

- **Stable Revenue Growth:** Maintaining a strong and consistent annual revenue over the last three years with an average growth rate of 10% per year.
- **Robust Capital Resources:** Supported by substantial charter capital and total assets, ensuring the capability for large-scale production expansion and significant inventory reserves.
- **Financial Transparency:** Fully audited financial reports are maintained and available upon request to verify our strong payment reputation and financial capacity.
- **High Operational Efficiency:** Consistently achieving healthy annual net profits with a stable growth trend; a gross profit margin of 28–32% reflects our superior manufacturing and business efficiency.

III. QUALITY MANAGEMENT SYSTEM

Our entire manufacturing system operates under a unified and rigorous quality management framework, fully compliant with ISO 13485 and ISO 9001 standards.



For every product in our portfolio, we maintain a comprehensive set of regulatory documentation:



- Trademark Registration Certificates:
- Product Testing Reports & Certificates of Conformity
- Certificates of Free Sale (CFS)
- Product Health Claim Registrations
- Technical Dossiers of Standards
- CO & CQ (Certificate of Origin & Certificate of Quality)

Notably, Tanaphar's Syringes and Needles have been granted the CE Marking.

IV. Global Export Capacity




- Tanaphar has successfully positioned itself as a reliable international supplier, bridging the gap between high-quality Vietnamese manufacturing and global healthcare needs.
- Global Footprint: We have successfully exported our medical consumables to over 30 countries worldwide.
- Market Leader in Asia: Tanaphar maintains a dominant market share across Southeast Asia and the broader Asian region, where our products are trusted for their consistent quality and competitive pricing.



PRODUCT PORFOLIO



No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
1		LA HAN QUA TANA Cough Lozenges	<p>Indications: Supports lung health, reduces phlegm, warms the throat; relieves cough due to cold, wind, dry cough, productive cough, weather changes, or persistent cough; reduces throat irritation.</p> <p>Ingredients per lozenge:</p> <ul style="list-style-type: none"> ❖ Luo Han Guo extract: 200mg; ❖ Terminalia chebula extract: 150mg; ❖ Ivy leaf extract: 100mg; ❖ Peppermint oil: 12mg; ❖ Plectranthus amboinicus oil: 9mg; ❖ Eucalyptus oil: 7mg. <p>Natural extracts, suitable for adults and children.</p>	3 Years	<ul style="list-style-type: none"> ❖ GMP ❖ Trademark Reg ❖ Testing Report ❖ Product Publication
2		LA HAN QUA TANA (Sugar-Free)	<p>Indications: Supports lung health, reduces phlegm, warms the throat; relieves cough due to cold, wind, dry cough, productive cough, weather changes, or persistent cough; reduces throat irritation.</p> <p>Ingredients per lozenge:</p> <ul style="list-style-type: none"> ❖ Luo Han Guo extract: 200mg; ❖ Terminalia chebula extract: 150mg; ❖ Ivy leaf extract: 100mg; ❖ Peppermint oil: 12mg; ❖ Plectranthus amboinicus oil: 9mg; ❖ Eucalyptus oil: 7mg. <p>Formulated for diabetics and people on restricted diets.</p>	3 Years	<ul style="list-style-type: none"> ❖ GMP ❖ Trademark Reg ❖ Testing Report ❖ Product Publication



No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
3		<p>EUCATANA Cough Lozenges</p>	<p>Indications: Supports lung health, soothes the throat, reduces phlegm, and warms the throat. Relieves cough due to cold, wind, dry cough, productive cough, or weather changes. Reduces throat pain and hoarseness.</p> <p>Ingredients per lozenge:</p> <ul style="list-style-type: none"> ❖ Peppermint oil: 9mg; ❖ Cinnamon oil: 1mg; ❖ Eucalyptus oil: 0.7mg; ❖ Ginger oil: 0.7mg; ❖ Plectranthus amboinicus oil: 0.5mg. <p>Natural extracts, suitable for adults and children.</p>	3 Years	<ul style="list-style-type: none"> ❖ GMP ❖ Trademark Reg ❖ Testing Report ❖ Product Publication
4		<p>STREPTANA Cough Lozenges</p>	<p>Indications: Supports lung health, soothes the throat, reduces phlegm, and warms the throat. Relieves cough due to cold, wind, dry cough, productive cough, or weather changes. Reduces throat pain and hoarseness.</p> <p>Ingredients per lozenge:</p> <ul style="list-style-type: none"> ❖ Peppermint oil: 7mg; ❖ Cinnamon oil: 1mg; ❖ Eucalyptus oil: 0.6mg; ❖ Ginger oil: 0.6mg. <p>Natural extracts, suitable for adults and children.</p>	3 Years	<ul style="list-style-type: none"> ❖ GMP ❖ Trademark Reg ❖ Testing Report ❖ Product Publication



No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
5		<p>UGOTANA Medical Adhesive Bandage</p>	<p>Usage: Protects small wounds, scratches, punctures, and minor skin tears.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Tape: Elastic Viscose and Polyamide fabric. ❖ Pad: Highly absorbent cotton with non-stick PE mesh. ❖ Adhesive: Solvent-free Zinc Oxide, Acrylic. 	3 Years	<ul style="list-style-type: none"> ❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report ❖ Standard Dossier
6		<p>UGOTANA Medical Adhesive Bandage</p>	<p>Usage: Protects small wounds, scratches, punctures, and minor skin tears.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Tape: Elastic Viscose and Polyamide fabric. ❖ Pad: Highly absorbent cotton with non-stick PE mesh. ❖ Adhesive: Solvent-free Zinc Oxide, Acrylic. 	3 Years	<ul style="list-style-type: none"> ❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report ❖ Standard Dossier
7		<p>UGOTANA Transparent Adhesive Bandage</p>	<p>Usage: Protects small wounds, scratches, punctures, and minor skin tears.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Tape: Polyethylene (PE). ❖ Pad: Highly absorbent cotton with non-stick PE mesh. ❖ Adhesive: Acrylic. 	3 Years	<ul style="list-style-type: none"> ❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report ❖ Standard Dossier




No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
8		UGOTANA Non-woven Adhesive Bandage	Usage: Protects small wounds, scratches, punctures, and minor skin tears. Components: ❖ Tape: Non-woven fabric. ❖ Pad: Highly absorbent cotton with non-stick PE mesh. ❖ Adhesive: Acrylic	3 Years	❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report ❖ Standard Dossier
9		UGOTANA Kids Adhesive Bandage	Usage: Protects small wounds, scratches, punctures, and minor skin tears. Components: ❖ Tape: Polyethylene (PE) with fun patterns. ❖ Pad: Highly absorbent cotton with non-stick PE mesh. ❖ Adhesive: Acrylic.	3 Years	❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report ❖ Standard Dossier
10		SUREMED Medical Silk Tape (Roll)	Usage: Breathable, easy to tear. Used for surgical dressings, fixing catheters, infusion needles, and medical devices. Material: ❖ Silk fabric base with solvent-free ❖ Zinc Oxide adhesive. Hypoallergenic, suitable for sensitive skin.	3 Years	❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report ❖ Standard Dossier




No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
11		<p>UGOTANA Medical Silk Tape (Roll)</p>	<p>Usage: High breathability, easy to tear. Used for surgical dressings, fixing catheters, infusion needles, and other medical devices.</p> <p>Material:</p> <ul style="list-style-type: none"> ❖ Silk fabric base with solvent-free ❖ Zinc Oxide adhesive. <p>Hypoallergenic, suitable for sensitive skin.</p>	3 Years	<ul style="list-style-type: none"> ❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report ❖ Standard Dossier
12		<p>UGOTANA Sterile Gauze Dressing</p>	<p>Usage: Covers and protects post-surgical wounds, scratches, skin tears, minor burns, and skin grafts. Can be used as a secondary dressing.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Soft breathable non-woven tape; ❖ Absorbent pad with non-stick PE mesh; ❖ Acrylic adhesive. <p>Sterile individual packaging.</p>	5 Years	<ul style="list-style-type: none"> ❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report ❖ Standard Dossier



No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
13		<p align="center">UGOTANA Waterproof Sterile Gauze Dressing</p>	<p>Indications: Used to protect post-surgical wounds, scratches, skin tears, minor burns, and skin grafts. Can be used as secondary dressing for other types of bandages.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Tape: Polyurethane (PU) film, thin, transparent, waterproof. ❖ Gauze: Highly absorbent cotton, surface covered with non-stick Polyethylene mesh. ❖ Adhesive: Acrylic, hypoallergenic, latex-free, leaves no residue upon removal. <p><i>The product is transparent for easy observation of fluid absorption; sterilized and individually packed to ensure sterility.</i></p>	5 Years	<ul style="list-style-type: none"> ❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report ❖ Standard Dossier
14		<p align="center">UGOTANA Waterproof Sterile Dressing</p>	<p>Indications: Used to protect shallow, dry wounds and thin skin areas frequently subject to friction. Also helps fix gauze, IV needles, tubes, central and peripheral venous catheters, while creating a barrier to protect the wound from external agents and bacterial penetration.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Tape: Polyurethane (PU) film, thin, transparent, waterproof. ❖ Adhesive: Acrylic, hypoallergenic, latex-free, leaves no residue upon removal. <p><i>The product is transparent for easy observation of the dressing area; sterilized and individually packed to ensure sterility.</i></p>	5 Years	<ul style="list-style-type: none"> ❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report ❖ Standard Dossier




No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
15		<p align="center">UGOTANA Waterproof Sterile Umbilical Gauze Dressing</p>	<p>Indications: Helps shield and keep the umbilicus clean and dry, preventing dirt and bacteria from entering, and limiting friction with diapers and clothes. This helps the umbilicus dry quickly, fall off naturally, and reduces the risk of dangerous infections.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ - Tape: Thin Polyurethane (PU) film, transparent, waterproof. ❖ - Gauze: Highly absorbent cotton, surface covered with non-stick Polyethylene mesh. ❖ - Adhesive: Acrylic, hypoallergenic, latex-free, leaves no residue upon removal. <p><i>The product is sterilized and individually packed, ensuring sterility and safety during use.</i></p>	5 Years	<ul style="list-style-type: none"> ❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report <p>Standard Dossier</p>
16		<p align="center">UGOTANA Sterile C-Section Gauze Dressing</p>	<p>Indications: Helps shield and protect the surgical site from dirt and bacteria, limiting painful collisions. Simultaneously, the dressing supports fast drying, good scarring, and reduces the risk of infection.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Tape: Non-woven fabric, soft, breathable, highly elastic. ❖ Gauze: Highly absorbent cotton, surface covered with non-stick Polyethylene mesh. ❖ Adhesive: Acrylic, hypoallergenic, latex-free, leaves no residue upon removal. <p><i>The product is sterilized and individually packed, ensuring sterility and safety during use.</i></p>	5 Years	<ul style="list-style-type: none"> ❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report ❖ Standard Dossier




No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
15		<p>UGOTANA Waterproof Sterile C-Section Gauze Dressing (10cm x 20cm)</p>	<p>Indications: Helps shield and protect the surgical site from dirt and bacteria, limiting painful collisions. Supports fast drying and reduces the risk of infection.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Tape: Polyurethane (PU) film, transparent, waterproof. ❖ Gauze: Highly absorbent cotton, surface covered with non-stick Polyethylene mesh. ❖ Adhesive: Acrylic, hypoallergenic. <p><i>The product is transparent for easy observation of the dressing area; sterilized and packed to ensure sterility.</i></p>	5 Years	<ul style="list-style-type: none"> ❖ ISO 13485:2016 ❖ Trademark Reg ❖ Testing Report ❖ Standard Dossier
16		<p>SOFTTANA Baby Cotton Swabs (Ear & Nose)</p>	<p>Indications: Helps gently clean the ear and nose areas for babies; effectively absorbs mucus and dirt without causing scratches or irritation.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Cotton head: 100% natural cotton fiber. ❖ Stem: Flexible medical plastic, high elasticity. <p><i>Safe for use on infants.</i></p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report



No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
17		<p>Bông vệ sinh tai mũi Người lớn SOFTTANA</p>	<p>Indications: Cleans dirt, mucus, nasal discharge, and earwax inside and around the ears and nose, aiding personal hygiene and preventing ear-nose issues. Also used for cleaning small details, makeup, or applying medicine.</p> <ul style="list-style-type: none"> ❖ Components: ❖ Cotton head: 100% natural cotton fiber. ❖ Stem: Flexible medical plastic, high elasticity 	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report
18		<p>TANAPHAR Sterile Surgical Gauze (8 layers / 12 layers)</p>	<p>Indications: Used for hemostasis, dressing wounds and surgical incisions; highly safe and convenient.</p> <p>Ingredients:</p> <ul style="list-style-type: none"> ❖ Absorbent gauze made of 100% cotton, neutral pH, free from starch or Dextrin. ❖ Water-soluble substances: not exceeding 0.5%. 	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
19		<p>TANAPHAR Alcohol Prep Pads</p>	<p>Indications: Used for cleaning and disinfecting the skin prior to injection.</p> <p>Ingredients:</p> <ul style="list-style-type: none"> ❖ 70% Isopropyl Alcohol. ❖ Absorbent cotton. 	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier

No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
20		<p>Chi nha khoa DENTANA</p>	<p>Indications: Helps eliminate plaque and food debris in interdental spaces where conventional toothbrushes cannot reach.</p> <p>Ingredients:</p> <ul style="list-style-type: none"> ❖ PTFE fiber. ❖ Flavoring. <p><i>Daily use ensures oral hygiene, strong teeth, and prevents bad breath.</i></p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
21		<p>ORALTANA Dental Floss</p>	<p>Indications: Helps eliminate plaque and food debris in interdental spaces where conventional toothbrushes cannot reach.</p> <p>Ingredients:</p> <ul style="list-style-type: none"> ❖ PTFE fiber. ❖ Flavoring. <p><i>Daily use ensures oral hygiene, strong teeth, and prevents bad breath.</i></p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
22		<p>ORALTANA Dental Floss Picks</p>	<p>Indications: A perfect combination of dental floss and toothpicks, specifically designed to eliminate plaque and food debris in interdental spaces where conventional toothbrushes cannot reach.</p> <p>Ingredients:</p> <ul style="list-style-type: none"> ❖ Floss: Nylon fiber. ❖ Handle: High-standard virgin plastic. <p><i>Daily use ensures oral hygiene, strong teeth, and prevents bad breath</i></p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier




No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
23		<p>ORALTANA 5-Star Dental Floss Picks (Grape / Mint / Peach flavors)</p>	<p>Indications: A perfect combination of dental floss and toothpicks, specifically designed to eliminate plaque and food debris in interdental spaces where conventional toothbrushes cannot reach.</p> <p>Ingredients:</p> <ul style="list-style-type: none"> ❖ Floss: Nylon fiber. ❖ Handle: High-standard virgin plastic. ❖ Flavoring: Grape, Mint, Peach. <p><i>Daily use ensures oral hygiene, strong teeth, and prevents bad breath.</i></p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
24		<p>QUICKTANA HCG Pregnancy Test Strip (5mm)</p>	<p>Indications: Helps detect early pregnancy 7 to 10 days after intercourse by identifying the presence of Human Chorionic Gonadotropin (HCG) hormone in urine.</p>	30 Months	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
25		<p>AMESTICK HCG Pregnancy Test Strip</p>	<p>Indications: Helps detect early pregnancy 7 to 10 days after intercourse by identifying the presence of Human Chorionic Gonadotropin (HCG) hormone in urine.</p>	30 Months	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier



No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
26		QUICKSEVEN HCG Pregnancy Test Strip	Indications: Helps detect early pregnancy 7 to 10 days after intercourse by identifying the presence of Human Chorionic Gonadotropin (HCG) hormone in urine.	30 Months	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
27		QUICKTANA HCG Pregnancy Test Pen (Midstream)	Indications: Helps detect early pregnancy 7 to 10 days after intercourse by identifying the presence of Human Chorionic Gonadotropin (HCG) hormone in urine.	30 Months	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
28		LH-TANA LH Ovulation Test Strip	Indications: Qualitative LH test strip helps determine ovulation day by detecting the increase in Luteinizing Hormone (LH) concentration in urine. When LH reaches its peak concentration, ovulation will occur after 24 hours and before 48 hours.	30 Months	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier




No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
29		<p>OVUTANA LH Ovulation Test Strip 5mm</p>	<p>Indications: Qualitative LH test strip helps determine ovulation day by detecting the increase in Luteinizing Hormone (LH) concentration in urine. When LH reaches its peak concentration, ovulation will occur after 24 hours and before 48 hours.</p>	30 Months	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
30		<p>HANOMED Nasal Aspirator</p>	<p>Indications: A baby's nasal cavity is easily blocked due to phlegm, mucus, milk, or food (due to choking). In such cases, the baby will have difficulty breathing and cry continuously, which can sometimes lead to danger. The nasal aspirator helps quickly clean and clear the nasal cavity and airway without causing pain to the baby.</p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
31		<p>SAKOMI Nasal Aspirator</p>	<p>Indications: A baby's nasal cavity is easily blocked due to phlegm, mucus, milk, or food (due to choking). In such cases, the baby will have difficulty breathing and cry continuously, which can sometimes lead to danger. The nasal aspirator helps quickly clean and clear the nasal cavity and airway without causing pain to the baby.</p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier

No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
32		<p>SOFTTANA Silicone Soft-tip Nasal Aspirator</p>	<p>Indications: A baby's nasal cavity is easily blocked due to phlegm, mucus, milk, or food (due to choking). In such cases, the baby will have difficulty breathing and cry continuously, which can sometimes lead to danger. The nasal aspirator helps quickly clean and clear the nasal cavity and airway without causing pain to the baby.</p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
33		<p>TANAKO Silicone Soft-tip Nasal Aspirator</p>	<p>Indications: A baby's nasal cavity is easily blocked due to phlegm, mucus, milk, or food (due to choking). In such cases, the baby will have difficulty breathing and cry continuously, which can sometimes lead to danger. The nasal aspirator helps quickly clean and clear the nasal cavity and airway without causing pain to the baby.</p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier


No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
34		<p>TANAKO Plus Silicone Soft-tip Nasal Aspirator</p>	<p>Indications: A baby's nasal cavity is easily blocked due to phlegm, mucus, milk, or food (due to choking). In such cases, the baby will have difficulty breathing and cry continuously, which can sometimes lead to danger. The nasal aspirator helps quickly clean and clear the nasal cavity and airway without causing pain to the baby.</p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
35		<p>TANA Silicone Nasal Irrigation Nozzle for Babies</p>	<p>Indications: The nasal irrigation nozzle helps prevent and limit ear, nose, and throat diseases. The product is made from soft, flexible virgin resin, fits snugly in the nostrils, providing a comfortable feeling and not causing damage to the mucosa. By reducing water pressure from the syringe, the nozzle helps clean and clear the nasal cavity quickly, gently, and painlessly.</p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
36		<p>TANAPHAR Sterile Pediatric Medical Tongue Cleaners</p>	<p>Indications: Helps clean the oral cavity, tongue, gums, and teeth, preventing oral diseases such as: tongue thrush, oral candidiasis, bad breath, gingivitis, tooth decay, etc.</p> <p>Components: Made from pure, ultra-fine, and soft fibers.</p> <p>Packaging: 01 product per individual pouch, sterilized to ensure absolute hygiene.</p>	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier




No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
37		<p>Oraltana Sterile Oral & Tongue Cleaners</p>	<p>Indications: Helps clean the oral cavity, tongue, gums, and teeth, preventing oral diseases such as: tongue thrush, oral candidiasis, bad breath, gingivitis, tooth decay, etc.</p> <p>Components: Made from pure, ultra-fine, and soft fibers.</p> <p>Packaging: 01 product per individual pouch, sterilized to ensure absolute hygiene.</p>	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
38		<p>Oraltana Sterile Oral & Tongue Cleaners - Wet</p>	<p>Indications: Helps clean the oral cavity, tongue, gums, and teeth, preventing oral diseases such as: tongue thrush, oral candidiasis, bad breath, gingivitis, tooth decay, etc.</p> <p>Components: Made from pure, ultra-fine, and soft fibers.</p> <p>Packaging: 01 product per individual pouch, sterilized to ensure absolute hygiene.</p>	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
39		<p>TANAKO Breast Pump</p>	<p>Indications: Supports breastfeeding mothers when the infant cannot latch directly, through its ability to express and store breast milk. The product is especially useful in cases of inverted, flat, or sore nipples; blocked milk ducts; oversupply of milk; or when the mother is away from the baby and milk cannot be expressed due to ineffective nursing.</p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier




No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
40		<p>TANAPHAR 3-Layer Medical Face Mask (Sterilized)</p>	<p>Indications: Helps prevent odors, smoke, fine dust, and bacteria from entering; supports the prevention of respiratory diseases.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Air filtration layer: Waterproof, breathable, non-allergenic to the skin. ❖ Nose wire: Easy to shape, helping to fix the mask and create a seal, fully covering the face. ❖ Earloops: Sturdy, with good elasticity for a comfortable fit. <p><i>Individually packed and EO sterilized.</i></p>	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
41		<p>TANAPHAR 3-Layer Medical Face Mask</p>	<p>Indications: Helps prevent odors, smoke, fine dust, and bacteria from entering; supports the prevention of respiratory diseases.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Air filtration layer: Waterproof, breathable, non-allergenic to the skin. ❖ Nose wire: Easy to shape, helping to fix the mask and create a seal. ❖ Earloops: Sturdy and elastic. 	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier




No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
42		<p>TANAPHAR 4-Layer Medical Face Mask</p>	<p>Indications: Provides enhanced protection to prevent odors, smoke, fine dust, and bacteria; supports the prevention of respiratory diseases.</p> <p>Components: 4 layers including high-quality non-woven fabric and filtration layers. Waterproof, breathable, and skin-friendly.</p>	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
43		<p>TANAPHAR Activated Carbon Medical Face Mask (Sterilized)</p>	<p>Indications: Helps prevent odors, smoke, fine dust, and bacteria; supports the prevention of respiratory diseases.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Filtration & 100% pure activated carbon layers: Waterproof, breathable, and non-allergenic. ❖ Nose wire & Earloops: Designed for a secure and comfortable fit. <p><i>Individually packed and EO sterilized.</i></p>	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
44		<p>TANAPHAR Activated Carbon Medical Face Mask</p>	<p>Indications: Helps prevent odors, smoke, fine dust, and bacteria; supports the prevention of respiratory diseases.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Filtration & 100% pure activated carbon layers: Waterproof, breathable, and non-allergenic. ❖ Nose wire & Earloops: Designed for a secure and comfortable fit. 	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier




No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
45		<p>COOLTANA Fever Cooling Patch</p>	<p>Indications: Provides natural cooling through heat absorption and diffusion from the body to the outside; supports temperature reduction, fever relief, and seizure prevention; soothes symptoms such as headache, toothache, muscle pain, stress, or heatstroke.</p> <p>Ingredients:</p> <ul style="list-style-type: none"> ❖ Aluminum glycinate ❖ Menthol ❖ Glycerin ❖ Eucalyptol ❖ Sodium Polyacrylate ❖ Water. <p><i>Provides continuous cooling effect for up to 10 hours; safe for both children and adults.</i></p>	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
46		<p>COOLING SHEET APPLE Fever Cooling Patch</p>	<p>Indications: Provides natural cooling through heat absorption and diffusion from the body to the outside; supports temperature reduction, fever relief, and seizure prevention; soothes symptoms such as headache, toothache, muscle pain, stress, or heatstroke.</p> <p>Ingredients:</p> <ul style="list-style-type: none"> ❖ Aluminum glycinate ❖ Menthol ❖ Glycerin ❖ Eucalyptol ❖ Sodium Polyacrylate ❖ Water. <p><i>Provides continuous cooling effect for up to 10 hours; safe for both children and adults.</i></p>	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier




No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
47		<p>DOCTORCOOL Fever Cooling Patch</p>	<p>Indications: Provides natural cooling through heat absorption and diffusion from the body to the outside; supports temperature reduction, fever relief, and seizure prevention; soothes symptoms such as headache, toothache, muscle pain, stress, or heatstroke.</p> <p>Ingredients:</p> <ul style="list-style-type: none"> ❖ Aluminum glycinate ❖ Menthol ❖ Glycerin ❖ Eucalyptol ❖ Sodium Polyacrylate ❖ Water. <p><i>Provides continuous cooling effect for up to 10 hours; safe for both children and adults.</i></p>	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
48		<p>COOLTANA Color-Changing Fever Patch</p>	<p>Indications: Provides natural cooling to relieve fever, headaches, and muscle pain.</p> <p>Smart Technology: Utilizes superior color-changing technology to monitor body temperature easily. When the body temperature reaches 38.5°C, the patch changes color from Purple to Pink, allowing for fever identification without a thermometer.</p> <p><i>Continuous cooling for up to 10 hours. Contains no pharmacological active ingredients, ensuring safety.</i></p>	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier




No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
49		TANAPHAR Sterile Wooden Tongue Depressor	<p>Indications: Used for oral, tongue, and throat examinations. Sterilized and single-use only, ensuring patient safety against the transmission of dangerous oral diseases.</p> <p>Components: Made from natural wood with fine, uniform grain, no core, and no warping. The wood is steamed before production, then dried and polished as a finished product to achieve high rigidity and smoothness for safe use.</p>	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
50		SOFTOMY Disposable Colostomy Bag (Sizes Φ 45, Φ 60)	<p>Indications: Used to collect waste for individuals unable to perform normal bowel movements. Secure and hygienic.</p>	3 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier
51		SAKURA Digital Thermometer	<p>Indications: Used to measure body temperature via mouth, armpit, or rectum, providing fast and accurate results. Features a clear LCD screen for easy monitoring and an automatic shut-off mode after 9–10 minutes of inactivity to save battery and increase durability.</p>	Unlimited	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier



No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
52		<p>TANA KO Digital Thermometer</p>	<p>Indications: Used to measure body temperature via mouth, armpit, or rectum, providing fast and accurate results. Features a clear LCD screen for easy monitoring and an automatic shut-off mode after 9–10 minutes of inactivity to save battery and increase durability.</p>	Unlimited	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report Standard Dossier
53		<p>DURATANA Condoms</p>	<p>Indications: Helps prevent pregnancy while avoiding the transmission of HIV/AIDS and other sexually transmitted diseases. Manufactured from natural rubber according to ISO 13485:2016 standards. Each condom is electronically tested to ensure optimal safety and reliability.</p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate
54		<p>TANAPHAR Hypodermic Needles (Various sizes)</p>	<p>Indications: Used for injections, medicine withdrawal, and blood sampling. Made from stainless steel with a triple-beveled sharp tip and no burrs to reduce pain during use. The needle shaft is smooth, uniformly round, has good mechanical rigidity, no warping, no internal impurities, and is silicone-coated to increase lubricity and safety during injection.</p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration Certificate ❖ ISO 13485:2016 Certificate ❖ Product Testing Report ❖ Standard Dossier

No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
55		TANAPHAR Single-use Syringe 1ml/cc	<p>Indications: Used for medical injections or pumping</p> <ul style="list-style-type: none"> ❖ Barrel, Plunger, and Gasket: Manufactured from medical-grade virgin resin, ensuring transparency, rigidity, and safety. ❖ Needle: Made of stainless steel, triple-beveled sharp tip, burr-free; smooth and round shaft with good mechanical rigidity, no warping, and silicone-coated for easier and safer injection. 	5 Years	<ul style="list-style-type: none"> ❖ CE Marking (European Export Standard) ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report ❖ Standard Dossier
56		TANAPHAR Single-use Syringe 3ml/cc	<p>Indications: Used for medical injections or pumping</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Barrel, Plunger, and Gasket: Manufactured from medical-grade virgin resin, ensuring transparency, rigidity, and safety. ❖ Needle: Made of stainless steel, triple-beveled sharp tip, burr-free; smooth and round shaft with good mechanical rigidity, no warping, and silicone-coated for easier and safer injection. 	5 Years	<ul style="list-style-type: none"> ❖ CE Marking (European Export Standard) ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report ❖ Standard Dossier
57		TANAPHAR Single-use Syringe 5ml/cc	<p>Indications: Used for medical injections or pumping.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Barrel, Plunger, and Gasket: Manufactured from medical-grade virgin resin, ensuring transparency, rigidity, and safety. ❖ Needle: Made of stainless steel, triple-beveled sharp tip, burr-free; smooth and round shaft with good mechanical rigidity, no warping, and silicone-coated for easier and safer injection. 	5 Years	<ul style="list-style-type: none"> ❖ CE Marking (European Export Standard) ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report ❖ Standard Dossier

No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
58		TANAPHAR Single-use Syringe 10ml/cc	<p>Indications: Used for medical injections or pumping solutions.</p> <ul style="list-style-type: none"> ❖ Barrel, Plunger, and Gasket: Manufactured from medical-grade virgin resin, ensuring transparency, rigidity, and safety. ❖ Needle: Made of stainless steel, triple-beveled sharp tip, burr-free; smooth and round shaft with good mechanical rigidity, no warping, and silicone-coated for easier and safer injection. 	5 Years	<ul style="list-style-type: none"> ❖ CE Marking (European Export Standard) ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report ❖ Standard Dossier
59		TANAPHAR Single-use Syringe 20ml/cc	<p>Indications: Used for medical injections or pumping solutions.</p> <ul style="list-style-type: none"> ❖ Barrel, Plunger, and Gasket: Manufactured from medical-grade virgin resin, ensuring transparency, rigidity, and safety. ❖ Needle: Made of stainless steel, triple-beveled sharp tip, burr-free; smooth and round shaft with good mechanical rigidity, no warping, and silicone-coated for easier and safer injection. 	5 Years	<ul style="list-style-type: none"> ❖ CE Marking (European Export Standard) ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report ❖ Standard Dossier
60		TANAPHAR Single-use Syringe 50ml/cc	<p>Indications: Used for medical injections or pumping solutions.</p> <ul style="list-style-type: none"> ❖ Barrel, Plunger, and Gasket: Manufactured from medical-grade virgin resin, ensuring transparency, rigidity, and safety. ❖ Needle: Made of stainless steel, triple-beveled sharp tip, burr-free; smooth and round shaft with good mechanical rigidity, no warping, and silicone-coated for easier and safer injection. 	5 Years	<ul style="list-style-type: none"> ❖ CE Marking (European Export Standard) ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report ❖ Standard Dossier

No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
61		TANAPHAR Single-use Feeding Syringe 50ml/cc	Indications: Supports the delivery of milk, liquid food, or nutritional solutions into the patient's body, helping care for and nourish patients who have difficulty eating directly. Single-use design ensures hygiene, safety, and convenience. Barrel, plunger, and gasket are made from virgin resin with high transparency and rigidity.	5 Years	<ul style="list-style-type: none"> ❖ CE Marking (European Export Standard) ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report ❖ Standard Dossier
62		TANAPHAR / HANOMED Infusion Set	<p>Indications: Used for intravenous infusion in medical practice.</p> <p>Components:</p> <ul style="list-style-type: none"> ❖ Tubing: Soft, flexible, durable, with high elasticity; does not kink during storage and use. ❖ Needle: Made of stainless steel, triple-beveled sharp tip, burr-free; smooth shaft with good mechanical rigidity, silicone-coated for safer and smoother operation. ❖ 	5 Years	<ul style="list-style-type: none"> ❖ CE Marking (European Export Standard) ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report ❖ Standard Dossier
63		TANAPHAR / HANOMED Infusion Set with Butterfly Needle	<p>Indications: Used for medical infusion.</p> <p>Components: Soft, kink-resistant tubing. The needle is triple-beveled and silicone-coated. Specifically, the needle features "butterfly wings" to provide stable fixation on the skin, limiting movement during infusion and ensuring patient comfort and safety.</p>	5 Years	<ul style="list-style-type: none"> ❖ CE Marking (European Export Standard) ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report ❖ Standard Dossier

No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
64		TANAKO Infusion Set	Indications: Used for medical infusion. Components: Soft, flexible, durable tubing with high elasticity. The needle features "butterfly wings" for stable fixation on the skin, minimizing movement and providing comfort for the patient. Silicone-coated for smooth operation.	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration ❖ ISO 13485:2016
65		TANAKO Butterfly Needle Set	Indications: Used for injections and infusions. Made from stainless steel with a triple-beveled sharp tip. The "butterfly wings" ensure secure fixation on the skin, limiting movement during the procedure and ensuring patient safety and comfort. Silicone-coated for smooth entry.	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration ❖ ISO 13485:2016
66		TANAPHAR Suction Catheter	Indications: Helps clear secretions and open the airway, facilitating air exchange. Supports the collection of secretions for diagnostic purposes while preventing infection due to stagnant phlegm. When used for deep suctioning (lower respiratory tract), the tube also stimulates the cough reflex, helping the patient expectorate secretions more effectively.	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report ❖ Standard Dossier

No	Product Image	Product Name	Indications & Ingredients	Shelf life	Certifications
67		TANAPHAR Newborn Umbilical Cord Clamp (Sterile - Single Use)	<p>Indications: Used to clamp the umbilical cord after cutting, helping to secure and protect the newborn's navel safely. The product prevents the penetration of bacteria, dirt, and harmful external agents, reducing the risk of umbilical infection – a common cause of dangerous complications in infants.</p> <p>Components: Manufactured from medical-grade virgin resin. Features a smooth surface to prevent skin injury, ensuring durability and a secure grip while remaining gentle and safe for newborns.</p> <p><i>Sterilized, single-use, providing maximum hygiene and convenience in neonatal care.</i></p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report ❖ Standard Dossier
68		TANAPHAR Plastic Vaginal Speculum (Sterile - Single Use)	<p>Indications: Used in gynecological examinations to assist physicians in observing and diagnosing reproductive health conditions conveniently and safely.</p> <p>Components: Made from transparent medical-grade virgin resin with a smooth surface to avoid mucosal injury during use.</p> <p><i>The speculum is sterilized by EO (Ethylene Oxide) gas, ensuring absolute sterility and eliminating the risk of cross-infection. The single-use design ensures maximum hygiene and convenience for both physicians and patients.</i></p>	5 Years	<ul style="list-style-type: none"> ❖ Trademark Registration ❖ ISO 13485:2016 ❖ Product Testing Report ❖ Standard Dossier

(*) NOTE:

- ❖ Manufacturing Origin: All products are wholly manufactured at the Tanaphar Joint Stock Company facility.
- ❖ Customization (OEM/ODM): Product specifications and packaging can be customized according to specific customer requirements and market preferences.

CORE COMMITMENTS TO GLOBAL PARTNERS

With over 20 years of manufacturing expertise and a facility infrastructure that meets global benchmarks, Tanaphar guarantees a seamless and high-quality supply chain for our international distributors.

- **Full Process Autonomy:** We maintain 100% control over the end-to-end process—from raw material sourcing and manufacturing to sterilized packaging and QC. This autonomy ensures rapid lead times and strictly adheres to both domestic and international medical regulations.
- **Supply Stability & Competitive Edge:** Our high production capacity ensures an abundant and stable supply even during market surges, backed by a pricing strategy designed to optimize our partners' market competitiveness.
- **Strategic Marketing Support:** Tanaphar is dedicated to our partners' success. We provide comprehensive support for in-store displays, promotional campaigns, and localized marketing communications to drive sell-through and elevate the end-user experience.

Direct contact:

- Alex Dao (Ms.) |Supply Chain & Import-Export Director
- TANAPHAR JOINT STOCK COMPANY
- A: Thanh Oai Industrial Complex, Binh Minh Commune, Hanoi City, Vietnam
- M: (+84) 968 174 279
- W: <http://www.tanaphar.com>
- F: <https://www.facebook.com/Tanapharjscimage.jpeg>